

Round Snap-Lock Containers Vented

Our 16 oz clear round plastic snap-lock polypropylene(PP) containers are now vented to allow for optimal gas exchange and feature tamper-evident security, letting you and your customers know if a product has been compromised.

The included tamper-evident lids are easy to open and use.

Round Snap-Lock Containers Vented Pricing

25 pk
100 pk
250 pk
500 pk

Specifications:

- Capacity: 16 oz (473 ml)
- Material Group: Plastic
- Material Type: Polypropylene (PP)
- Color: Transparent
- Shape: Round
- Diameter: 3.6 in
- Height: 3.2 in
- Weight: 21.7 grams
- Cap: Included



Vented Plant TC Test Tubes

The Vented Plant Tissue Culture (TC) Test Tubes (80ml) are a product designed for modern plant tissue culture laboratories and hobbyists.

Caps and rack included. The test tubes come in packs of 21.

Vented Plant TC Test Tubes Pricing

21 pk

Specifications:

- Diameter (width) of the base: 1.25 in
- Height: 6 inches
- Included: 21 Vented test tubes, vented caps, plastic rack

